

**Hydrophilic Single-use Catheter** Men Boys Women Girls

Single sterile packages. Length without funnel. Available in packs of 30.

	Size*	no additional wrapping		with additional wrapping	
		CH	mm	Order no./REF	Order no./REF
<b>Functional Packaging</b> IQ-Cath 20	length 40 cm, not ready for use, functional packaging, sterile				
		10	3,3	IQ 2004.10	-
		12	4,0	IQ 2004.12	IQ 2014.12
		14	4,7	IQ 2004.14	IQ 2014.14
		16	5,3	IQ 2004.16	-
<b>Functional Packaging</b> IQ-Cath 21	length 40 cm, ready for use, sterile saline solution 0.9 % integrated, functional packaging, sterile				
		10	3,3	IQ 2104.10	IQ 2114.10
		12	4,0	IQ 2104.12	IQ 2114.12
		14	4,7	IQ 2104.14	IQ 2114.14
		16	5,3	IQ 2104.16	IQ 2114.16
<b>"Peel" Packaging</b> IQ-Cath 29	length 40 cm, not ready for use, "Peel" packaging, sterile				
		10	3,3	IQ 2904.10	-
		12	4,0	IQ 2904.12	-
		14	4,7	IQ 2904.14	-
		16	5,3	IQ 2904.16	-
<b>"Peel" Packaging</b> IQ-Cath 29	length 12 cm, not ready for use, "Peel" packaging, sterile				
		10	3,3	IQ 2902.10	-
		12	4,0	IQ 2902.12	-
		14	4,7	IQ 2902.14	-
		16	5,3	IQ 2902.14	-
<b>"Peel" Packaging</b> IQ-Cath 23	length 40 cm, <b>immediately ready for use</b> , sterile saline solution 0.9 % activated, "Peel" packaging, sterile				
		10	3,3	IQ 2304.10	-
		12	4,0	IQ 2304.12	-
		14	4,7	IQ 2304.14	-
		16	5,3	IQ 2304.16	-
<b>"Peel" Packaging</b> IQ-Cath 23	length 30 cm, <b>immediately ready for use</b> , sterile saline solution 0.9 % activated, "Peel" packaging, sterile				
		10	3,3	IQ 2303.10	-
		12	4,0	IQ 2303.12	-
		14	4,7	IQ 2303.14	-
		16	5,3	IQ 2303.14	-
<b>"Peel" Packaging</b> IQ-Cath 23	length 12 cm, <b>immediately ready for use</b> , sterile saline solution 0.9 % activated, "Peel" packaging, sterile				
		10	3,3	IQ 2302.10	-
		12	4,0	IQ 2302.12	-
		14	4,7	IQ 2302.14	-
		16	5,3	IQ 2302.14	-
<b>"Peel" Packaging</b> IQ-Cath 43	length 40 cm, <b>immediately ready for use</b> , sterile saline solution 0.9 % activated, Sleeve (insertion aid), "Peel" packaging, sterile				
		10	3,3	IQ 4304.10	-
		12	4,0	IQ 4304.12	-
		14	4,7	IQ 4304.14	-
		16	5,3	IQ 4304.16	-
<b>"Peel" Packaging</b> IQ-Cath 43	length 30 cm, <b>immediately ready for use</b> , sterile saline solution 0.9 % activated, Sleeve (insertion aid), "Peel" packaging, sterile				
		10	3,3	IQ 4303.10	-
		12	4,0	IQ 4303.12	-
		14	4,7	IQ 4303.14	-
		16	5,3	IQ 4303.14	-
<b>Urine Collection Bag</b> IQ-Cath 26	length 40 cm, ready for use, sterile saline solution 0.9 % integrated, urine collection bag 1 litre, sterile				
		12	4,0	IQ 2604.12	IQ 2614.12
		14	4,7	IQ 2604.14	IQ 2614.14
		16	5,3	IQ 2604.16	IQ 2614.16
		18	6,0	IQ 2604.18	-





\*The size of the catheter refers to its external diameter indicated in French gauge (CH) and millimetres. The funnel is coloured to code the respective size.

**Single-use Gel Catheter**



Single sterile packages. Length without funnel. Available in packs of 30.

Size*		no additional wrapping	with additional wrapping
CH	mm	Order no./REF	Order no./REF

<b>Functional Packaging</b>	<b>IQ-Cath 30</b>	length 40 cm, not ready for use, functional packaging, sterile		♂	○ 12 4,0	<b>IQ 3004.12</b>	-
					● 14 4,7	<b>IQ 3004.14</b>	-
					● 16 5,3	<b>IQ 3004.16</b>	-
<b>Functional Packaging</b>	<b>IQ-Cath 31</b>	length 40 cm, ready for use, gel integrated, functional packaging, sterile		♂	● 10 3,3	<b>IQ 3104.10</b>	-
					○ 12 4,0	<b>IQ 3104.12</b>	<b>IQ 3114.12</b>
					● 14 4,7	<b>IQ 3104.14</b>	<b>IQ 3114.14</b>
					● 16 5,3	<b>IQ 3104.16</b>	<b>IQ 3114.16</b>
<b>Functional Packaging</b>	<b>IQ-Cath 32</b>	length 40 cm, <b>immediately ready for use</b> , prelubricated with gel, functional packaging, sterile		♂	● 10 3,3	<b>IQ 3204.10</b>	-
					○ 12 4,0	<b>IQ 3204.12</b>	-
					● 14 4,7	<b>IQ 3204.14</b>	-
					● 16 5,3	<b>IQ 3204.16</b>	-
			● 18 6,0	<b>IQ 3204.18</b>	-		
<b>Urine Collection Bag</b>	<b>IQ-Cath 36</b>	length 40 cm, ready for use, gel integrated, urine collection bag 1 litre, sterile		♂	○ 12 4,0	<b>IQ 3604.12</b>	<b>IQ 3614.12</b>
					● 14 4,7	<b>IQ 3604.14</b>	<b>IQ 3614.14</b>
					● 16 5,3	<b>IQ 3604.16</b>	-

\*The size of the catheter refers to its external diameter indicated in French gauge (CH) and millimetres. The funnel is coloured to code the respective size.

